

A man and a woman are shown from the back, interacting with a large, white medical device. The man is shirtless and the woman is wearing a yellow t-shirt. The device has a control panel with a small screen and several buttons. The background is a plain, light-colored wall.

pausch
m e d i c a l

PAUSCH Medical Product Portfolio – Distinguishing Features and Product Attributes






*Made in
Germany*

Your Partner for High-Quality Medical Technology





For more than 80 years, PAUSCH Medical GmbH has produced equipment for the medical industry and in particular for X-ray diagnostics. With our highly dedicated and experienced development team, we deliver individualized product solutions for our customers, making PAUSCH Medical a long sought-after partner among the global leaders in medical device manufacturing.

Our product portfolio includes tabletop, column-, and ceiling-mounted systems in the X-ray sector, with numerous accessories. The devices are designed for operation with analogue systems as well as with digital detectors, and can be flexibly retrofitted as needed.

Universal X-Ray Units

	Cosmos2	Disy Med	Patienten Positioning Table (movable)
			
Product description	The universal X-ray system enables image recording of either standing or seated patients, and of recumbent patients with the optional patient positioning table. With multiple setting options, no additional X-ray system is necessary.	A single detector solution with unique flexibilities, and a multifunctional motorized rotatable column system. In contrast to columns with straight carrier arms, lateral images on recumbent patients are also possible using Disy Med.	Free wheeled patient table for easy and effortless positioning.
Dimensions	height: 244.6 cm // width: 226 cm // depth: 130 cm	height: 320 cm // width: 224.1 cm // depth: 142.4 cm	length: 200 / 185 cm // width: 65 cm // height: 71 cm
Weight	approx. 550 kg	approx. 600 kg (incl. add-on components)	57 kg
Patient bearing load			210 kg
Floating tabletop			
Fixed height			✓
Height adjustable			
Vertical	✓	✓	
Tiltable	✓	✓	
X-Ray Generator installable			
Digital	✓	✓	
Analogue	✓	✓	
Fluoroscopy			
Tomography			
Accessories	✓	✓	✓
CE	✓	✓	✓
UL	✓	✓	on request

Bucky Tables

	Elevator 3000	Covator	CS3000	CS2000 Tomo	Combi-Elevator Tomo
					
Product description	Height adjustable bucky table with floating table top. Usable together with our ceiling mount tube stand and our bucky wall stand.	A comprehensive solution for small X-ray rooms including a height-adjustable Bucky with an integrated column system and wall mount.	The bucky table with integrated column system as a compact solution in connection with a bucky wall stand.	Bucky table with integrated column system and linear tomography.	Height adjustable bucky table with integrated column system and linear tomography.
Dimensions (length x depth x height)	220 x 104 x 55 - 90 cm	235 x 138.4 x 227.5 cm	220 x 139.7 x 227.5 cm	229.5 x 151.5 x 233.5 cm	229.5 x 151.5 x 233.5 cm
Weight	approx. 250 kg	approx. 600 kg	approx. 400 kg	approx. 420 kg	approx. 540 kg
Patient bearing load	320 kg	320 kg	450 kg	228 kg	228 kg
Floating tabletop	✓	✓	✓	✓	✓
Fixed height			✓	✓	
Height adjustable	✓	✓			✓
Vertical					
Tiltable					
X-Ray Generator installable		✓	✓		
Digital	✓	✓	✓	✓	✓
Analogue	✓	✓	✓	✓	✓
Fluoroscopy					
Tomography				✓	✓
Accessories	✓	✓	✓	✓	✓
CE	✓	✓	✓	✓	✓
UL	✓	in process	on request	✓	✓

3D Ceiling Systems

Ceiling Mount Tube Stand TL3000



Product description

Motorized ceiling system with integrated tracking functions. Variable ceiling rails, automatic format size sensing and stitching modes as further options. Together with our Elevator and our WS3000 it is possible to configure all options and site combinations fully free.

Dimensions height: 94.5 - 246.5 cm // depth: 127 cm

Weight 320 - 360 kg

Patient bearing load

Floating tabletop

Fixed height

Height adjustable ✓

Vertical ✓

Tiltable ✓

X-Ray Generator installable

Digital ✓

Analogue ✓

Fluoroscopy

Tomography

Accessories ✓

CE in process

UL in process


Bucky Wall Stands

BS2000 V ⁽¹⁾

Optional: BS2000 VK ⁽²⁾

WS3000 V ⁽¹⁾

Optional: WS3000 VK ⁽²⁾



A standard wall-mounted unit for versatile usage.

Space-saving design with the same high stability. The WS3000 can be installed as free-standing without external braces, and is equipped with electromagnetic brakes.

height: 221.9 cm // width: 65.1 cm // depth: 42.5 cm

141 kg

height: 224.5 cm // width: 65.1 cm // depth: 36.7 cm

141 kg

✓

Optional: BS2000 VK ⁽²⁾

✓

Optional: WS3000 VK ⁽²⁾

✓

✓

✓

✓

✓

✓

✓

✓

✓

in process


⁽¹⁾ Vertical adjustable ⁽²⁾ Vertical adjustable and tiltable

Urology

Uromat 3000



Uromat DFD (Flat Panel)



Product description

A urological examination and operating table, which is equipped with cutting-edge digital X-ray technology. It is suitable for urological examinations and diagnostic work, as well as for minimally invasive urological surgery. This includes transdermal procedures, such as endoscopies, kidney and urinary stone operations, and procedures involving the prostate, bladder, and ureter.

The Uromat DFD provides convenient access for the patient as well as optimum and ergonomic operating conditions for the user. By using dynamic digital X-Ray detectors of up to 43 x 43 cm (17" x 17") you're capturing the entire urinary (kidney, ureter, bladder) in one single shot.

Dimensions footprint: 370 x 210 cm // height: 260 cm

footprint: 370 x 210 cm // height: 260 cm

Weight approx. 900 kg

approx. 900 kg

Patient bearing load 205 kg

250 kg

Floating tabletop ✓

✓

Fixed height

Height adjustable ✓

✓

Vertical ✓

✓

Tiltable ✓

✓

X-Ray Generator installable

Digital ✓

✓

Analogue ✓

Fluoroscopy ✓

✓

Tomography

Accessories ✓

✓


CE ✓

✓

UL on request

✓



A dark blue world map serves as the background for the top half of the page. A central white circle in Europe is connected by dashed white lines to six other white circles located in North America, South America, Africa, Asia, and Australia. The text is positioned on the left side of the map.

PAUSCH Medical GmbH
Graf-Zeppelin-Strasse 1
91056 Erlangen
Germany

Phone +49 (0)9131/99 92 - 0
Fax +49 (0)9131/99 24 - 22
eMail info@pausch.de

© Copyright 2015
PAUSCH Medical GmbH, Erlangen, Germany
Subject to change.

www.pauschmedical.com

A company of **MAVIG GROUP**